

40

Forty Years of
Swiss Innovation



Spot on ...

LED Contra Angle 20:1

Art.-No. 5052

20:1 LED contra angle with hexagon chucking system

Precision is your daily business, we have the suitable instrument with light to achieve the best results on the working spot.



We bring light into dark

The contra angle with a ratio of 20:1 is a standard implantology instrument with Intra coupling for torque driven motor systems. The built in LED-Light gives the clear view you need in daily business.

- Dynamo LED light system
- Dismantleable in 3 parts for easy cleaning
- Natural white light for realistic color rendering
- Internal and external cooling system
- Light independent from control unit due to integrated dynamo



Technical Data, 20:1 LED contra angle

Type:	20:1 LED contra angle
Speed:	up to 40,000 rpm
Illumination:	5,000 – 31,000 lux
Instruments shank diameter:	2.35 mm
Chuck system:	Push button, hexagon
Length:	110 mm
Weight:	100 g

Distribution partner:



Nouvag AG · St.Gallerstr. 23-25 · CH-9403 Goldach
Phone +41 (0)71 846 66 00 · Fax +41 (0)71 846 66 70
info@nouvag.com · www.nouvag.com